


HPLC-06

LabVIEW routine to automatically acquire chromatograms. Acquisition is triggered from the HPLC/UHPLC autosampler.

We created this virtual instrument using LabVIEW for the acquisition of signals (chromatograms) from our HPLC and UHPLC systems. It uses USB devices from National Instruments. The program can be triggered either manually or automatically (by the autosampler). It and it stores the chromatograms as IGOR data that can be opened with our HPLC analysis (see its description on this web).



This program is continuously running in our lab for amine analyses and can be freely shareware upon request.